

SecurSeal® Replaces IVA® Seals

USP <797> Compliant

Zero Paper — Zero Cardboard
No Shedable Material — Compare Pricing & Save

SecurSeal® IV Seals

Clear Poly Release Liner

- Perforated for Easy Use
- No Paper
- No Shedable Material

Plastic Core

- No Cardboard
- No Shedable Material

Clear all Plastic Dispenser Box

- No Cardboard
- Wipe Down Compatible
- No Shedable Material

Available From:

**Cardinal
McKesson
Williams Medical
AmerisourceBergen
Morris & Dickson
H.D. Smith**



CALL TOLL FREE
800.527.4616
or visit us at www.winfieldlabs.com



See wholesaler cross reference chart on the other side.

IVA® Seals have paper liner and cardboard core that generate particulate matter, as shown on the black tape.




IVA® is a trademark of Covidien AG or an affiliate

SecurSeal® IV Seals Wholesaler Cross Reference Chart

SERIES 10

Bag Seal
Baxter Bag
1000/Roll




If you currently use CP3011

| Color/SecurSeal | Cardinal | AmerisourceBergen | McKesson |
|-----------------|----------|-------------------|----------|
| Silver 10-100 | 3518156 | 993232 | 1181684 |
| Blue 10-200 | 3518164 | 993257 | 1172261 |
| Red 10-300 | 3518172 | 993269 | 1179316 |
| Green 10-400 | 3518180 | 993271 | 1175553 |
| Yellow 10-500 | 3518198 | 993283 | 1192640 |
| Orange 10-600 | 3518206 | 993290 | 1177997 |
| Tan 10-700 | 3518214 | 993307 | 1185495 |

Fits Hospira VisIV® Container and other large port bags

SERIES 30

Large
Bag Seal
1000/Roll




If you currently use CP 3013

| Color/SecurSeal | Cardinal | AmerisourceBergen | McKesson |
|-----------------|----------|-------------------|----------|
| Silver 30-100 | 3563947 | 993795 | 1126523 |
| Blue 30-200 | 3563970 | 993802 | 1203314 |
| Red 30-300 | 3563996 | 993814 | 1212976 |
| Green 30-400 | 3564002 | 993826 | 1205558 |
| Yellow 30-500 | 3564028 | 993838 | 1142918 |
| Orange 30-600 | 3564036 | 993845 | 1210137 |
| Tan 30-700 | 3564044 | 993853 | 1140326 |

SERIES 28

28mm
Bottle Tops
1000/Roll




If you currently use CP 3001

| Color/SecurSeal | Cardinal | AmerisourceBergen | McKesson |
|-----------------|----------|-------------------|----------|
| Silver 28-100 | 3562675 | 993720 | 1185818 |
| Blue 28-200 | 3562683 | 993737 | 1175579 |
| Red 28-300 | 3562733 | 993749 | 1182385 |
| Green 28-400 | 3562758 | 993752 | 1178219 |
| Yellow 28-500 | 3563780 | 993764 | 1198522 |
| Orange 28-600 | 3562782 | 993776 | 1179324 |
| Tan 28-700 | 3562816 | 993788 | 1194083 |

SERIES 55

Hospira Bag
Side Injection
Port Bag 1000/Roll




If you currently use CP 3004

| Color/SecurSeal | Cardinal | AmerisourceBergen | McKesson |
|-----------------|----------|-------------------|----------|
| Silver 55-100 | 3562840 | 993434 | 1212786 |
| Blue 55-200 | 3562857 | 993446 | 1198407 |
| Red 55-300 | 3562907 | 993473 | 1208826 |
| Green 55-400 | 3562915 | 993480 | 1203199 |
| Yellow 55-500 | 3563764 | 993509 | 1140201 |
| Orange 55-600 | 3562956 | 993511 | 1205483 |
| Tan 55-700 | 3562980 | 993547 | 1124148 |

SERIES 20

20mm
Vial Tops
1000/Roll




Also fits
other large
port bags

If you currently use CP 3002

| Color/SecurSeal | Cardinal | AmerisourceBergen | McKesson |
|-----------------|----------|-------------------|----------|
| Silver 20-100 | 3563012 | 993555 | 1168608 |
| Blue 20-200 | 3563038 | 993562 | 1156439 |
| Red 20-300 | 3563046 | 993574 | 1164318 |
| Green 20-400 | 3563905 | 993586 | 1160381 |
| Yellow 20-500 | 3563913 | 993612 | 1172691 |
| Orange 20-600 | 3563921 | 993624 | 1163294 |
| Tan 20-700 | 3563939 | 993636 | * |

SERIES 13

13mm
Vial Tops
1100/Roll



If you currently use CP 3003

| Color/SecurSeal | Cardinal | AmerisourceBergen | McKesson |
|-----------------|----------|-------------------|----------|
| Silver 13-100 | 3563806 | 993648 | 1153139 |
| Blue 13-200 | 3563830 | 993651 | 1142884 |
| Red 13-300 | 3563848 | 993670 | 1151984 |
| Green 13-400 | 3563855 | 993687 | 1144492 |
| Yellow 13-500 | 3563863 | 993699 | 1154400 |
| Orange 13-600 | 3563871 | 993701 | 1149111 |

If your wholesaler shows "not stocked" on any item, Winfield will drop ship items until your distribution center begins stocking that item number. Free shipping - No Minimum required!



Call Toll Free **800.527.4616** Or Fax Toll Free **877.455.1665**
Visit Us at www.winfieldlabs.com